

Specifications of Pneumatic Compressor Device

S.N	Parameters	Technical Specifications	Remarks
1	Usage	For the prevention of Deep Vein Thrombosis in a bed ridden Patient	
2	Power Supply	220 V, 50 Hz	
3	Features	Dual Cuff Design	
		Inflation time maximum 12 seconds & deflation time maximum 18 seconds.	
		LED display for pressure on specific limb which glows when cuff on that limb is inflated.	
		Knobs for each cuff for setting & adjusting the pressure.	
		Audio & visual safety alarms for leak	
4	Accessories	In the event of high pressure cuff should deflate automatically.	
		Cuffs should have skin friendly cotton lining.	
		Colour coded cuff sizes (Medium, large, XL, XXL, XXXL) (each size cuff should be supplied)	
5	Terms and Conditions	Free installation and demo at Pokhara MH and Itahari MH	
		Warranty : Minimum 2 years/Provided by OEM whichever is greater	
		Manufacture Date must be less than 2 years at time of delivery.	
		Should be backed with reliable service.	
		Standard : ISO/CE Equivalent	

Handwritten notes:
 1. *Handwritten signature*
 2. *Handwritten signature*
 3. *Handwritten signature*
 4. *Handwritten signature*
 5. *Handwritten signature*

Handwritten notes:
 1. *Handwritten signature*
 2. *Handwritten signature*
 3. *Handwritten signature*
 4. *Handwritten signature*
 5. *Handwritten signature*

Handwritten notes:
 1. *Handwritten signature*
 2. *Handwritten signature*
 3. *Handwritten signature*
 4. *Handwritten signature*
 5. *Handwritten signature*